

## Product Details

With suppliers in N/A

# TabTimer Med-e-lert Pill Dispenser (PR01241)

<https://ilcaustralia.org.au/products/1241>

This battery operated circular unit alerts the user when medication is due and dispenses the dosage. The alert is either an audible alarm or flashing light. It contains 28 compartments and enables up to six alarms per day. The unit is lockable and only the right medication is exposed for dispensing.

## General Details

### Caution: Medication Management

Always check with a pharmacist before storing medications in a dispenser or storage container as some medications require special storage conditions and some medications cannot be kept in the same container as others.

### Safety and Usage Advice

Refer to the supplier and manufacturer manual for maintenance instructions and safety warnings.

### Extra product information

Med-O-Wheel Smart can be programmed for one to four dose times daily as required, and it reminds the patient to take medication on time with audio and light alarms. Once the patient pushes the transparent part of the device cover, it rotates the cassette inside and makes the correct dose available, while waiting for the next dose time. If a dose has been missed despite all the reminders, the dispenser will skip that dose and make the next correct dose available, when its time comes. The missed dose remains inside the cassette to allow that medication to be monitored and re-used.

### Method of Use

1. Press button 1 for 3 seconds - the time starts blinking.
2. Set hours with button 2 and minutes with button 3. Accept by pressing button 1.
3. C:28/14 indicates the number of compartments. Change by pressing button 3. Accept by pressing button 1.
4. L:15-0 indicates audio volume. Change by pressing button 3.



*PR01241 TabTimer Med-e-lert Pill Dispenser*

Accept by pressing button 1.

5. Set dose times, 1 - 4 times / day. Set hours with button 2 and minutes with button 3. NB! If medication is to be taken twice a day, for example, set dose time 3 to --:-- ( no time ). Accept by pressing button 1.

6. Press button1 for 3 seconds.

7. Press button 3 for 3 seconds. The display will read SI-HO. Button 1 moves the cassette one notch, button 2 returns the cassette to the initial position. Press button 3 for 3 seconds. The display then shows the time and number of compartments ( 14 or 28 ).

When the lid is closed, press on the display to check the next dose time.

### **Manufacturer/Country**

China

## Models / Specifications

### **Model Med-e-Lert Pill Dispenser: One Model Only**

- **Depth Overall**
- **Width Overall**
- **Number of compartments**
- **Diameter** 193 mm
- **Depth** 52 mm
- **Weight Overall** 48 g

## Suppliers (0)

## Agents (0)

## Importers (0)

## Manufacturers (0)

## Maintenance

Refer to the supplier and manufacturer manual for maintenance instructions and safety warnings.

## Standards

Some assistive technology needs to meet Australian or other standards. Standards may relate

to materials, manufacturing, and installation. Products that meet Australian or international standards will have written certification. To find out if a product meets Australian Standards ask the supplier to show you the certificate. For more information about standards see [http://www.ilcaustralia.org.au/Using\\_Assistive\\_Technology/standards\\_regulation](http://www.ilcaustralia.org.au/Using_Assistive_Technology/standards_regulation)

## Disclaimer

Independent Living Centres Australia (ILCA) does not design or make products. ILCA cannot guarantee that all information given is correct. Prices shown are a rough guide. Some products may incur GST. Check with suppliers for up-to-date prices. The display of products at an ILCA centre, or details given in print, is only for your information. It is not a recommendation. You should always get professional advice when choosing products. The person choosing the product is responsible for their choice. ILCA takes no responsibility for people choosing equipment or services which will cause (or be alleged to have caused) loss or injury. You must follow the manufacturer's or service provider's instructions when using equipment or services.

Copyright ILCA 2011. All rights reserved