

Product Details

With suppliers in N/A

Rifton Pacer Gait Trainer (PR05531)

<https://ilcaustralia.org.au/products/5531>

A rigid height adjustable gait trainer made of welded steel and aluminium for use with children and adults, which can be used indoor and outdoor. It has a variety of accessories that can be positioned to support the individual according to their individual needs. The frame is lightweight and foldable. It is available in five sizes. The frame size is determined by the users elbow height.

General Details

Options

- Arm prompt set
- Chest prompt
- Hip positioner
- Thigh prompt
- Ankle prompt

Extra product information

The frame can be used for gait training in anterior (forward facing) and posterior (backward facing). The frame height is adjusted in 25mm increments using a trigger.

There are four large castors for indoor and outdoor use. The castors can be individually locked to prevent swiveling, to prevent rolling backward, to control speed and to brake the wheels.

Accessories can be attached to the frame in the required position using clamps which are tightened using knobs and adjusted using press buttons on the post.

The frames are available in champagne, red, blue, purple, lime and raspberry.

===

The frame consists of welded steel and aluminium tubing with a baked on finish, adjustable trunk support and handle bars. The sliding components consist of polymer coated steel and aluminium. Straps are nylon and polypropylene. Upholstery is vinyl covered urethane foam.

The frame is able to be positioned in two ways. The user facing the open end of the frame (posterior positioning), or the closed end of the frame (anterior positioning).



*PR05531 Rifton Pacer Gait
Trainer*

The castors have five functions - gradual brake/drag, brake lock, swivel, swivel lock and one way ratchet control enabling control of speed, direction and level of manoeuvrability. Each castor can be locked individually to prevent swivelling.

The truck support components are padded and can be adjusted independently. The prompts are attached to the frame by clamps and do not require tools. Prompts are purchased separately.

Manufacturer/Country

Rifton, USA

Websites

User Manual - <https://www.rifton.com/~media/files/rifton/product-manuals/pacers-small-mini-na49.pdf>

Models / Specifications

Model Pacer Mini: User elbow height 390-520mm

- **Weight** 3.5 kg
- **Load Capacity (supplier stated)** 22.6 kg
- **Width Overall** 520 mm
- **Length Overall** 570 mm
- **Height of Armrest from Ground** 390 mm - 520 mm

Model Pacer Small: User elbow height 470-700mm

- **Weight** 6 kg
- **Length Overall** 830 mm
- **Width Overall** 570 mm
- **Load Capacity (supplier stated)** 34 kg
- **Height of Armrest from Ground** 470 mm - 700 mm

Model Pacer Large: User elbow height 790-1170mm

- **Width Overall** 700 mm
- **Length Overall** 970 mm
- **Load Capacity (supplier stated)** 91 kg
- **Weight** 13 kg
- **Height of Armrest from Ground** 790 mm - 1170 mm

Model Pacer Medium: User elbow height 610-860mm

- **Width Overall** 660 mm
- **Load Capacity (supplier stated)** 68 kg
- **Length Overall** 810 mm
- **Weight** 7 kg
- **Height of Armrest from Ground** 610 mm - 860 mm

Model Pacer X-Large: User elbow height 860-1270mm

- **Load Capacity (supplier stated)** 160 kg
- **Weight** 20 kg
- **Length Overall** 1140 mm
- **Width Overall** 810 mm
- **Height of Armrest from Ground** 860 mm - 1270 mm

Suppliers (o)

Agents (o)

Importers (o)

Manufacturers (o)

Images



PR05531 Rifton Pacer Gait Trainer

Maintenance

Refer to the supplier and manufacturer manual for maintenance instructions and safety warnings.

Standards

Some assistive technology needs to meet Australian or other standards. Standards may relate to materials, manufacturing, and installation. Products that meet Australian or international standards will have written certification. To find out if a product meets Australian Standards ask the supplier to show you the certificate. For more information about standards see http://www.ilcaustralia.org.au/Using_Assistive_Technology/standards_regulation

Disclaimer

Independent Living Centres Australia (ILCA) does not design or make products. ILCA cannot guarantee that all information given is correct. Prices shown are a rough guide. Some products may incur GST. Check with suppliers for up-to-date prices. The display of products at an ILCA centre, or details given in print, is only for your information. It is not a recommendation. You should always get professional advice when choosing products. The person choosing the product is responsible for their choice. ILCA takes no responsibility for people choosing equipment or services which will cause (or be alleged to have caused) loss or injury. You must follow the manufacturer's or service provider's instructions when using equipment or services.

